HD gas purification for polarized HDice target production at Jefferson Lab

C. S. Whisnant, A. D'Angelo, L. Colaneri, J. Devilbiss, T. Kageya, D. A. Loving, M. Lowry, A. Rizzo, A. M. Sandorfi, C. Schaerf, E. Speiser, J. D. Storey, C. M. Wallace, X. Wei, and I. Zonta

Overview

- Commercial HD Gas contains 0.5-2.0% H₂ and D₂
- To produce targets, the H₂ and D₂ must be reduced
 Require small, carefully controlled concentrations
- Isotopic separation is done using distillation
 - Operate at low temperature to liquefy and obtain good boiling point (partial pressure) differences
 - Apply heat to boil the liquid
 - Use a packed column for good efficiency
- Gas monitoring
 - Residual Gas Analyzer
 - Gas chromatography
 - Raman Scattering

Isotopic Separation



H₂, HD, and D₂ have significantly different vapor pressures at low temperature.

This allows separation using distillation techniques.





Mechanical refrigerator

- Two-stage system
 - 1st stage: ~35 K
 - 2nd stage

- Temperature control & boilup power supply
 - Typical boil-up power ~ 2.3W
- 18K → P ~ 260 mbar
- Extraction rate = 1 mole/day

Packed Distillation Column

- Heli-pak
- Improves efficiency
 - Nucl. Inst. Meth. 664 (2012)
 347
 - <u>http://www.wilmad-</u>
 <u>labglass.com/Products/</u>
 <u>LG-6730-104/</u>



Windows in vacuum container and thermal shield permit viewing the boiling liquid hydrogen

... has often enabled system debugging



Gas Handling

- Manifold with 2-mole stainless tanks
- Gas moved between tanks using LHe and a cold trap
- Long term storage at low pressure in Al tanks







Low Temperature Gas Chromatography



Chromatography Column

- Off-the-shelf PLOT 5A capillary column
 - ValcoPLOT[®] VP-Molesieve GC Capillary Column
 - Length=30 m, ID=0.53 mm, Film Thickness=50µm
 - <u>http://www.sisweb.com/gc/vici/vpmol.htm</u> (CFS-X-03053-500)
- Column operated between T≈115K 130K
 - Immersed in a roughly equal mixture of LN₂ and isopentane
 - ΔT during a run is 1-2K
 - Run-to-run variation is larger as LN₂ boils away
 - As T increases: separation and retention time decreases

Monitoring with Gas Chromatography



Raman spectroscopy: laser scattering from rotational states

Raman spectroscopy does not dissociate HD molecules and it is used to provide measurements of H_2 and D_2 relative concentrations in HD gas at the 10-5 level.



Molecular rotational energy levels are given by:

$$E_R = \hbar^2 \frac{J^2}{2I} = hcb_0 J(J+1)$$

where:

$$b_0 = \frac{n}{8\pi^2 I}$$

L

and I is the molecule moment of inertia.

Distance among levels is higher for lighter molecules.

Raman scattering occurs when Laser light is absorbed and re-emitted by a molecule while the rotational state changes: the scattered light is shifted from the Laser energy by:

$$\Delta E = hcb_0(J+2)(J+3) - hcb_0J(J+1) = hcb_0(4J+6)$$

Selection rules $\Delta J=\pm 2$ connect ortho- to ortho- and para- to para- molecular states.

Raman set-up

- The gas to be analyzed is used to fill the cell.
- The Laser light is focused at the center of the cell and reflected back by a mirror to increase the useful incoming power.
- The scattered light is reflected by a mirror and expanded in a parallel beam to enter the monochromator slit.



The scattered light spectrum is detected by a cooled CCD and recorded by a PC.









Raman spectrum analysis

- The relative content of H₂ and D₂ contaminants in HD gas is determined fitting all the peaks in the frequency range between 100 and 1000 cm-1 higher than the original Laser frequency (stokes lines).
- Peaks intensities are determined from the Gaussian fit parameters
- By studying the dependence upon the temperature and upon the index J of each rotational level it is possible to determine the relative number of molecules in the sample for the different chemical species.
- The minimum relative amount that may be detected is lowered by increasing the Laser power and by reducing the background fluctuations

we use a 15 Watts Laser and we integrate the spectrum over 1 hour to average over noise fluctuations

Comparing Raman and GC

- Raw HD before purification:
 - Raman
 - H2/HD = (7.0 \pm 0.6) X 10⁻³
 - D2/HD = (8.0 ± 0.07) X 10⁻³
 - -GC
 - H2: (7.0 ± 0.14) X 10⁻³
 - D2: (11.0 ± 0.22) X 10⁻³



H₂ and D₂ peaks were not visible in the chromatogram

Summary

- Low-temperature distillation permits separation of H_2 , HD, and D_2 to a few parts in 10⁻⁴ or better
- Measurement techniques available for monitoring purity (in increasing sensitivity and complexity)
 - Residual gas analysis
 - Gas chromatography
 - Raman scattering
- Gas production underway

The Need for Purity

Creating a frozen spin target with significant relaxation time requires that H_2 and D_2 concentrations be small and carefully controlled.

Longitudinal relaxation time T₁H

$$\frac{1}{T_1^H} \propto \left(c_1^H\right)^{-2.1} \left(c_1^D\right)^{-1.6}$$





Gas is stored at < 2 ATM in specially treated Al cylinders to minimize this effect.